TempoTIM
TempoTIM Plus

EN INSTRUCTIONS FOR USE

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INSTRUCTIONS FOR USE

Pedal exerciser, foldable

Pedal exerciser with calories counter, foldable

Product code: **TempoTIM / TempoTIM Plus**


1. Characteristics

The lower limb pedal exerciser is made of a powder-coated steel frame. It is equipped with resistance mechanism, four anti-slip pads and foot straps. It has a rotary knob with a wide range of resistance adjustment, making it suitable for various patient groups. It can be used in hospitals, hospices and at home as a rehabilitation and personal recreation device.

Elements of the exerciser:

- a. Foot protection strap
- b. Adjustment knob
- c. Locking pin
- d. Anti-slip pads
- e. Counter (model **TempoTIM Plus**)

2. Application

The pedal exerciser is designed for patients with contracture or muscular atrophy caused by long-term immobilisation due to a plaster bandage or surgery, as well as for patients who need to restore motor functions and re-learn the correct movement pattern due to damage to the neurological system. It is possible to adjust the exerciser to the state of limb fitness and to the progress in rehabilitation owing to resistance adjustment.

3. Assembly and adjustment

To assemble pedal exerciser, extend both legs until the locking pins click.

- Assemble adjustment knob into the mounting opening in the upper part of the frame.
- To fold the exerciser, push in the locking pin sequentially on both sides of the frame, then fold the exerciser legs inwards.
- To remove the counter or replace the battery in the counter (**TempoTIM Plus** model), completely unscrew the adjustment knob, then remove the counter and/or replace the battery.
- To replace the battery, release the battery cover latch by pressing the button on the bottom of the counter.

ATTENTION: All data will be lost after battery change.

Resistance adjustment

In order to increase fitness levels and improve muscle work and muscle strength gains, the training intensity and/or duration of exercise should be gradually increased.

- To increase the resistance, turn the adjustment knob clockwise.
- To reduce the resistance, turn the adjustment knob anticlockwise.
- When performing low impact exercises, the adjustment knob should be completely unscrewed, whereas when performing resistance exercises, the adjustment knob should be tightened to the desired load.

Usage

- Place the pedal exerciser on the floor and make sure the anti-slip pads are firmly in place to prevent the pedal exerciser from moving.
- Exercises should be carried out in a seated position, the pedal exerciser should be positioned exactly in front of the user so that the hip joint is bent at an angle of no less than 90°.

- Place feet on the pedals and adjust the strap to foot circumference, this will prevent the foot from slipping off the pedal during exercise.
- Set the resistance using the adjustment knob.
- If the exerciser has a counter (**TempoTIM Plus** model), select the corresponding function.

Counter function:

<p>CNT - Counter Activity counter</p>	
<p>SCAN - Function scan Display of exercise parameters</p>	
<p>T.CNT – Counter (total) Total number of rotation of all exercises</p>	
<p>TIME - Time Duration of exercise</p>	
<p>CAL - Calories Counter of calories burnt</p>	
<p>Auto Power ON - Automatic start Automatic start of the counter when activity is detected</p>	
<p>Auto Power OFF - Automatic stop Automatic stop of the meter after 5 min. of inactivity</p>	
<p>Pause TIMER - Automatic pause Automatic stop of the counter after 3 sec. of inactivity (display of "STOP" in the upper left-hand corner of the display)</p>	
<p>RESTART - Erase data Press and hold the red button for min. 2 sec. To erase T.CNT, remove the battery</p>	

4. Technical data

Base dimensions	39,5 cm x 47 cm
Base height from the ground	17 cm
Weight: TempoTIM / TempoTIM Plus	2,25 kg 2,30 kg
Maximum capacity	120 kg
Battery	1 x 1.5V (AAA)
Frame colour	black

5. Notes

Timago International Group does not bear responsibility of improper use of the pedal exerciser, failure to comply safety regulations and misuse.

6. Maintenance

Carry out regular maintenance of the technical condition of the pedal exerciser:

- It is recommended to clean the pedal exerciser regularly with a soft cloth and soapy water. After cleaning, dry the pedal exerciser thoroughly with soft cloth.
- Do not use strong polishing agents for cleaning as it may cause damage to the protective coatings.
- Remove impurities from moving elements of the exerciser on an ongoing basis.

7. Storage and transport conditions

Device should not be exposed to direct sunlight, water or moisture.

8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof purchase (receipt or invoice) should be retained.

9. Labels

	Reference numer		Manufacturer
	LOT numer		Manufacture date
	Serial numer		Medical device
	Note		Please read the manual
	The manufacturer checked for compliance with essential medical devices requirements.		

The entire range is available to you from our distributors.

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