



HK 44  
Cosy

beurer  
wellbeing




## HK 44 Cosy Heat pad - not a hot water bottle



Soft and cosy heat pad



Machine washable up to 30 °C



3 illuminated temperature  
settings



Soft and cosy heat pad in the shape of a hot water bottle

Breathable, cosy and skin-friendly

3 illuminated temperature settings

Automatic switch off after approx. 90 minutes

Rapid heat up

BSS (Beurer Safety System)

Electronic temperature control

Detachable control unit

Machine washable up to 30 °C

OEKO-TEX ®

The textiles used for this device meet the stringent human ecological requirements of OEKO-TEX ® Standard 100, as verified by Hohenstein Research Institute.

Mains operation: 220-240 V; 50-60 Hz; 100 W

Product dimensions (LxW): approx. 50 x 30 cm

Product weight: approx. 340 g

3 year guarantee

Sales unit: 5 / Shipping carton: -

EAN no.: 4211125 27505 1

Item no.: 275.05

