

Invacare Bed Cradle Statio L871

Data Sheet



Bed Cradle Statio L871 is one of Invacare's bed accessories designed to meet users' needs adequately. When you are confined to bed and no part of your body must be subject to any pressure, the bed-cradle supports the weight of the sheets and blankets. Bed cradle can be placed on the bed end to alleviate pressure from the duvet on the legs/ feet.

The Invacare Bed Cradle Statio L871 can help if your skin is sensitive. It helps with air circulation to your feet keep your skin dry. If you have burns, open skin sores, or infections, the Invacare Bed Cradle Statio L871 can also be quite useful. Bed Cradle Statio L871 is also ideal if you have certain injuries or conditions, such as paraplegia or a pressure injury. Also, if you lie in bed for long periods, Bed Cradle can also be useful. So, you can see the Bed Cradle Statio L871 can be applied in many areas.

Technical Data:

- Safe working Load 25kg

Material:

- Steel tube

Manufacturer:

Invacare Portugal, Lda Rua Estrada
Velha 949 4465-784 Leça do Balio
Portugal



Yes, you can: