

99766774 Normasol

Normal saline for use in topical irrigation, wound cleansing, irrigation of mucous membranes and rinsing of eyes

Product details

Size : 25ml
 Descriptive feature : Saline, Sodium chloride
 Sterile : Sterile

Images



Delivered items

99766774-02

Sales released in: Bahrain, Gibraltar, Hong Kong, Ireland, Malta, Oman, Saudi Arabia, South Africa, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Country of origin: United Kingdom of Great Britain and Northern Ireland

Shelf life: 1.5 years

Sterilization method: Steam

Production Responsibility: Mölnlycke Health Care Ltd, Tubiton House, Medlock Street, Oldham, OL1 3HS, United Kingdom

Packing information: First packaging layer is a plastic sachet. Once opened the sterile barrier cannot be closed again. Second layer is a thermoformed pouch. Third layer is a corrugated board transport box.

Is suitable for Tray: No

Packing level	Quantity	GS1 Code / UDI-DI	Width x Length x Height	Vol	Weight gross / net
Consumer pack	1	7332551878603			
Transport box	250	7323190024803	295x395x185 mm	21.6 dm3	7.6 / 6.6 kg
Pallet	16000	7323190024797	800x1200x1624 mm		

Material

Animal tissues : No

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00. The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2022 Mölnlycke Health Care AB. All rights reserved.



Human blood derivatives : No
 Natural rubber latex : No
 Medicinal substances : No
 Phthalates : No
 Polyvinyl chloride : No

Product Composition Saline Solution Products

Product Component	Composition
Saline solution	Purified water, Sodium chloride

Technical

Classifications

Regulation type(s)	MDR Class IIb	Locally Regulated	Locally Regulated	Locally Regulated
CE Certificate Number :	MDR 722028			
Intended Purpose :	For use in topical irrigation, wound cleansing, irrigation of mucous membranes and rinsing of eyes	For use in topical irrigation, wound cleansing, irrigation of mucous membranes and rinsing of eyes	For use in topical irrigation, wound cleansing, irrigation of mucous membranes and rinsing of eyes	For use in topical irrigation, wound cleansing, irrigation of mucous membranes and rinsing of eyes
MDR Classification Rule :	4			
Conformity Annexes :	IX			
Measuring Function :	No			
Notified body medical devices/PPE :	BSi 2797	Saudi Authority	Local Authority MEA	Local Authority APAC
510(k) clearance number :	3110756			
Regulatory Released :	Gibraltar, Ireland, Malta, United Kingdom of Great Britain and Northern Ireland	Saudi Arabia	Bahrain, Oman, South Africa	Hong Kong

Find out more at www.molnlycke.com



Classifications

Regulation type(s)	Locally Regulated
CE Certificate Number :	
Intended Purpose :	For use in topical irrigation, wound cleansing, irrigation of mucous membranes and rinsing of eyes
MDR Classification Rule :	
Conformity Annexes :	
Measuring Function :	
Notified body medical devices/PPE :	Local Authority United Arab Emirates
510(k) clearance number :	
Regulatory Released :	United Arab Emirates

Applied standards : The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 14001

Removable label

No

GMDN Code (Global Medical Device Nomenclature)

64871 Sterile wound irrigation sodium chloride solution, single-use

UNSPSC

42312309 Wound irrigation systems

Commodity Code

3307900000 Sodium Chloride Solutions

CE Responsibility

Molnlycke Health Care AB, Gamlestadsvägen 3c, Box 13080, SE-402 52 Goteborg Sweden

Basic UDIDI

733243000000000062JW

Find out more at www.molnlycke.com

