

TurnAid Complete System

AS EASY AS IT GETS!

Item no: 1000

TurnAid is a medical device aid developed to mount on a care bed. It combines the two well-known basic functions bedrail and electrically driven turning system.



TurnAid is developed to optimize the handling related to the care of, moving and handling of bedridden patients, with severe physical disabilities in hospitals, care homes, institutions or in their own home.

TurnAid is also suitable for persons in risk of getting pressure ulcers, as the construction of the system allows moving of the patient in the sheet, which reduces the friction that occurs between the patient and the sheet.

TurnAid handles persons with a weight of up to 200 kg, and does so with maximum consideration to both patient and staff.

TurnAid is easily mounted on most care beds, as it comes with a built-in universal fitting, which allows for length and width adjustment. This enables TurnAid to comply with the Medical Directive's strict demands in this area.



Specifications:

Total weight of the system	46 kg.
Length of the system	190 cm
Width of the system	10 cm
Height of the system	57 cm
User weight minimum/max	40/200 kg
Maximum work load	240 kg
Noise level	>70db (A)
Connection voltage	230 V / 50 Hz
Inspection/service	Yearly
IP class	IPX0= Not suitable for outdoor use and not flush tight.
Classification rule	Class I according to enclosure IX in 93/42/EEC
Applied standards	BS EN 60601-1 :2016 (3ed) BS EN 60601-1-2:2015 (EMC) –Parts BS EN 60601-2-52:2010

TurnAid complete system includes:

- 2 Bedrails with integrated electric turning system (right and left)
- 1 Control box with power cable, motor cables, bracket and hand control
- 1 Turning sheet
- 1 Glide sheet

Accessories:

- Turning sheets
- Glide sheets
- Extra hand control
- Foot control
- Transport and mounting trolley
- Transport bags

The system is secured against misuse when both hand- and foot control are connected. If both are activated at the same time, none of them will operate.